











Consumables Specifications



| Items | Dimensions & Standards (if applicable) | Photos | Applications |
|--|--|--|--|
| <p>3M Absorbent Pad (6 pieces)</p> | <p>19" x 16"</p> |  | <p>Absorb chemical spills (2.6L per piece).</p> <p>Spill type: Acids/Alkalis/Alcohols, Chemicals, Flammable Liquids, Solvents - Hydrocarbon Chlorinated</p> |
| <p>Bio-hazard Tag (100 pieces Yellow Identification Tag)</p> | <p>4.7" x 2.4"</p> |  | <p>To label each bag of clinical waste with the yellow bio-hazard tag for identification and enhancing traceability.</p> <p>Click to see instructions on ordering and using.</p> |



| Items | Dimensions & Standards (if applicable) | Photos | Applications |
|---------------------------------------|--|---|--|
| Broken Glass Box (with 1 plastic bag) | 14" x 13" x 19" |  | For broken glass disposal, the box should be disposed of with the broken glass to protect the waste collectors from cut injury. |
| L-Red Clinical Waste Bag (50 pieces) | 24" X 42" High Density Polyethylene (HDPE) 0.75mm |  | For disposal of clinical waste except Group 3: Human and Animal Tissues. Note: Sharp box can be put inside a red bag before disposal. |


| Items | Dimensions & Standards (if applicable) | Photos | Applications |
|---|--|---|--|
| S-Red Clinical Waste Bag (50 pieces) | 24" X 24" High Density Polyethylene (HDPE) 0.75mm |  | For disposal of clinical waste except Group 3: Human and Animal Tissues. Note: Sharp box can be put inside a red bag before disposal. |
| L-Yellow Clinical Waste Bag (50 pieces) | 24" X 42" |  | For disposal of Group 3 Clinical Waste: Human and Animal Tissues |

| Items | Dimensions & Standards (if applicable) | Photos | Applications |
|------------|--|---|---|
| Cryogloves | <p>6" x 15"</p> <p><u>Standard</u> EU standards for CE Cat. III: EN511, EN 388, EN 420</p> |  | <p>Dispensing/transferring cryogenic liquids, removing samples from cryogenic liquids, any open cryogenic containers with the chance of exposure to liquids or splashes</p> |
| Cryo-Apron | <p>24" X 42"</p> <p><u>Standard</u> EU standards for CE Cat. III: EN 511, EN 388, EN ISO 13688</p> |  | <p>Dispensing/transferring cryogenic liquids, handling of cryogenic samples into or removing from storage vessels</p> |

| Items | Dimensions & Standards (if applicable) | Photos | Applications |
|--------------------|--|---|---|
| Face Shield | 8" x 15.5" Standard : ANSI Z87.1 |  | Protect the wearer's face from hazards such as flying objects, chemical splashes, liquid nitrogen splashes. |
| Plastic Spill Tray | 14.2*14.2*X14.2" |  | A catch tray for storage of inorganic chemical waste in a chemical waste cabinet |

| Items | Dimensions & Standards (if applicable) | Photos | Applications |
|--|---|---|--|
| Safety Glasses | <p>NA</p> <p><u>Standard</u> GB14866-2006. ANSI Z787+- 2010,, CE EN 166</p> |  | <p>General use for the protection of eyes from dust, flying debris and liquids splash.</p> <p><i>PS: For the frequent handling of large amounts of corrosives or chemical liquids, a Chemical Goggle should be used.</i></p> |
| <p>Spill Kit Refill, includes:</p> <p>Nitrile Rubber Glove x 1 Disposable Plastic Glove x 5 Disposable Plastic Bag x 3 3M Absorbent Pad x 6 Pan and Brush x 1 Safety Glasses x 1 Labels (corrosive, flammable, harmful, toxic)</p> | <p>NA</p> |  | <p>For replenishing the consumables of a spill response kit that will be used to handle minor chemical spillage.</p> |

| Items | Dimensions & Standards (if applicable) | Photo | Applications |
|---|--|---|---|
| <p>Spill Response Kit, includes:</p> <p>Plastic Box x 1 Nitrile Rubber Glove x 1 Disposable Plastic Glove x 5 Disposable Plastic Bag x 3 3M Absorbent Pad x 6 Pan and Brush x 1 Safety Glasses x 1 Labels (corrosive, flammable, harmful, toxic)</p> | <p>15" x 10.2" x 8.3"</p> |  | <p>A standard kit with essential consumables to handle minor chemical spillage.</p> <p>Chemical using laboratories should be equipped with this kit or an equivalent alternative.</p> |
| <p>Mercury Spill Response Kit</p> | <p>NA</p> |  | <p>For the handling of mercury spillage.</p> |

| Items | Dimensions & Standards (if applicable) | Photos | Applications |
|------------------|---|--|---|
| Steel Spill Tray | 15" x 15" x 5.7" |  | A catch tray for storage of organic chemical waste in a chemical waste cabinet. |